


# NomoLine-0™ Capnography Sampling Lines

効果的なカプノグラフィモニタリングを可能にする気道アダプタセットとカニューーラ



## NomoLine-O™ Sampling Line Options

Division	Product Description	Part Number	Length	Image
High Humidity or Long-term Airway Adapter Set	NomoLine-O HH Adult/Paediatric Airway Adapter Set, 2m, 25/Box	4449	2m	
	NomoLine-O HH Infant/Neonatal Airway Adapter Set, 2m, 25/Box	4450	2m	
High Humidity or Long-term Cannula with O2 tubing	NomoLine-O HH Adult Nasal CO2 Cannula with O2, 2m, 25/Box	4451	2m	
	NomoLine-O HH Paediatric Nasal CO2 Cannula with O2, 2m, 25/Box	4645	2m	
	NomoLine-O HH Adult Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4453	2m	
	NomoLine-O HH Paediatric Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4647	2m	
High Humidity or Long-term Cannula without O2 tubing	NomoLine-O HH Adult Nasal CO2 Cannula, 2m, 25/box	4643	2m	
	NomoLine-O HH Paediatric Nasal CO2 Cannula, 2m, 25/box	4644	2m	
	NomoLine-O HH Infant Nasal CO2 Cannula, 2m, 25/box	4664	2m	
	NomoLine-O HH Adult Nasal/Oral CO2 Cannula, 2m, 25/box	4452	2m	
	NomoLine-O HH Paediatric Nasal/Oral CO2 Cannula, 2m, 25/box	4646	2m	
	NomoLine-O HH Adult Single Nasal Prong CO2 Cannula, 2m, 25/box	4665	2m	

## NomoLine-O 仕様

長さ	2 m
サンプリング流量	50 ml/分
インライン細菌フィルター	BFE > 99.9%
チューブ内径	1 mm
酸素補助供給	最大 5 l/分
気道アダプタ死腔	
成人 / 小児	≤ 6 ml
幼児および幼児 / 新生児	≤ 0.7 ml
材質	天然ゴムラテックス不使用

本製品は医療従事者用です。適応、禁忌、警告、注意を含む詳細な使用方法については、取扱説明書をご覧ください。

カナダでは、一部品番の提供を行っていません。

製造販売業者  
マシモジャパン株式会社  
<http://www.masimo.co.jp>

本社  
〒169-0074 東京都新宿区北新宿2-21-1  
新宿フロントタワー 24階  
TEL: 03-3868-5201 FAX: 03-3868-5202

